

APPLICATION PROCEDURE

OLYMPUS AU2700/5400

α-AMYLASE

mod. IFCC

Liquid
2 Reagents

REF

Cont.

D94570	5 x 100 ml	4 x 100 ml 1 x 100 ml	Reagent 1 Reagent 2	4000 tests 125 µl/test
D94571	5 x 50 ml	4 x 50 ml 1 x 50 ml	Reagent 1 Reagent 2	2000 tests 125 µl/test
D96569	5 x 10 ml	4 x 10 ml 1 x 10 ml	Reagent 1 Reagent 2	400 tests 125 µl/test
D98485	5 x 3 ml	Calibrator	Diacal Auto	
D98481	12 x 5 ml	Control normal	Diacon N	
D98482	12 x 5 ml	Control abnormal	Diacon P	

1. Reagent and Sample Preparation:

The reagents are ready to use. For Urine: dilute urine 1 + 1 with distilled water.

2. Instrument settings:

a) Substrate Start:

Temperature: 37 °C
Test Name **AMY** Type **Serum** Operation **Yes**

Sample	Volume	<input type="text" value="2.0"/>	µl	Dilution	<input type="text" value="0"/>	µl
Reagents	R1 Volume	<input type="text" value="100"/>	µl	Dilution	<input type="text" value="0"/>	µl
	R2 Volume	<input type="text" value="25"/>	µl	Dilution	<input type="text" value="0"/>	µl
Wavelength	Pri	<input type="text" value="410"/>		Sec	<input type="text" value="660"/>	
Method	<input type="text" value="RATE"/>					
Reaction Slope	<input type="text" value="+"/> +					
Measuring point 1	First	<input type="text" value="20"/>		Last	<input type="text" value="27"/>	
Measuring point 2	First	<input type="text"/>		Last	<input type="text"/>	
Linearity	<input type="text" value="15"/> %					
No-Lag-Time	<input type="text"/>					
Pre-dilution Rate	<input type="text"/>					
Min OD	<input type="text" value="-0.100"/>			Max OD	<input type="text" value="1.200"/>	
Reagent OD Limit	First L	<input type="text" value="-0.100"/>		First H	<input type="text" value="1.40"/>	
	Last L	<input type="text" value="-0.100"/>		Last H	<input type="text" value="1.40"/>	
Dynamic Range	L	<input type="text" value="0"/>		H	<input type="text" value="1200"/>	
Correlation Factor	A	<input type="text" value="1"/>		B	<input type="text" value="0"/>	
On-board stability period	<input type="text" value="40"/> Days					
Value/Flag	#	Level L	#	Level H	#	#
Normal Ranges	Sex	Year	Month	Year	Month	
1	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
2	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
3	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
4	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
5	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
6	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
7	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	

APPLICATION PROCEDURE

Panic Value		L	H	Unit	
	#	<input type="text"/>	#	<input type="text"/>	U/L
Calibration Type		<input type="text" value="AB"/>	Formula	<input type="text" value="Y=AX+B"/>	Counts
				<input type="text" value="#"/>	
	Cal.No	OD	CONC	Factor OD-L	Factor OD-H
Point 1	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="*"/>	<input type="text"/>	<input type="text"/>
Point 2	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="*"/>	<input type="text"/>	<input type="text"/>
Point 3	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="*"/>	<input type="text"/>	<input type="text"/>
Point 4	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="*"/>	<input type="text"/>	<input type="text"/>
Point 5	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="*"/>	<input type="text"/>	<input type="text"/>
Point 6	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="*"/>	<input type="text"/>	<input type="text"/>
Point 7	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="*"/>	<input type="text"/>	<input type="text"/>
1-Point Cal. Point		<input type="text"/>			
MB Type Factor		<input type="text"/>		Calibration Stability Period	
				<input type="text"/>	

- #) Data entry by the user
- *) Enter calibration or standard value and position
- ***) The factor must be checked by a calibration serum.
- ***) For urine

NOTE:

These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing the tests.