

APPLICATION PROCEDURE

Α ΕΛΥΟΙΨΠΩΙΟΑ ΕΕΨΥÙÙΝΟΨΟΑΓ Δ

GAMMA GT

Szasz, stand. to IFCC

Liquid

2 Reagents

REF

Cont.

D95604 5 x 50 ml 4 x 50 ml Reagent 1 1000 tests
 1 x 50 ml Reagent 2 250 µl/test

D98485 5 x 3 ml Calibrator Diacal Auto
 D98481 12 x 5 ml Control normal Diacon N
 D98482 12 x 5 ml Control abnormal Diacon P

N. Reagent preparation:
 The reagents are ready to use.

2. Instrument settings:

Item Name	GGT
Units	[U/l]
Decimals	[0]
Calibration Type	[Linear]
1	[*]
2	[*]
3	[]
4	[]
5	[]
6	[]
Method Type	[RATE]
MainW. Length1	[405]
Subw. Length2	[660]
Method	[Szasz]
N-Range Male and Female	[#] – [#]
Reagent Kind	[Double]
Sample Volume	[8 µl]
Reagent 1 Volume	[200 µl]
Reagent 2 Volume	[50 µl]
Co-ordinate Y = [1.000] X + [0.000]	
Monitor Zero level point	[1]
Monitor Span	[3.000]
Third mix	[Off]
Water Blank	[On]

Read Main (Start-End)	[35] – [54]
Sub	[] – []
Absorbance Limit	[-3.000] – [3.000]
Linear Check	[90%]
Endpoint Limit	[0.000]
Factor: Blank Correction	[1.000]
Prozone	
Read	1 [] – []
	2 [] – []
	3 [] – []
Limit	2 [] – []
	3 [] – []
Auto Rerun SW	[Off]
Range (Result)	
Low-SW	[Off]
Serum	[0]
High-SW	[Off]
Serum	[Off]
Absorbance	
Low-SW	[Off]
High-SW	[Off]
Prozone-SW	[Off]

#) Data entry by the user
 *) Input Standard or Calibrator

NOTE:

These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing the tests