

FURUNO 180 / DAYTONA

CREATININE

"Mod. Jaffe"

Liquid., 2 Reagents

REF

Cont.

D95595	5 x 100 ml	4 x 100 ml R1 1 x 100 ml R2	2500 tests/kit 200µl /test
D98485	5 x 3 ml	Calibrator	Diacal Auto
D98481	12 x 5 ml	Control normal	Diacon N
D98482	12 x 5 ml	Control abnormal	Diacon P

1. Reagent preparation:

The reagent is ready to use.

2. Instrument settings:

RunMnr(F5)	ChmPrm(F6)	Calibration (F7)	QC (F8)	System Prms	Maintenance								
Chemistry Parameters Method 11 Name CREA		R1 Reagent Name CREA Volume 160 uL R2 <input type="radio"/> Disable <input checked="" type="radio"/> Enable	Reagent CRE2 Volume 40 uL Name <input checked="" type="radio"/> Disable <input type="radio"/> Enable Wash Enable	Reagent Name Reagent Type Wash									
Unit µmol/L	Assay Type Rate	Measuring Points 1 <input checked="" type="radio"/> Disable 2	Start - End - Start 15 Decimal Points End 21 Normal Range	Dilution <input checked="" type="radio"/> Disable <input type="radio"/> Enable 0 53 - 115	Reagent Name Reagent Type Wash								
Wave Length Prim. 510 Sec. <input type="checkbox"/>	Sampling Volume 8 uL Dilution <input checked="" type="radio"/> Disable <input type="radio"/> Enable	Technical Range (Conc.) 0 - 1330 (mAbs/10) - 50000	RPT Wash (R1) (R2) <input checked="" type="radio"/> SysWater <input type="radio"/> Wash Sol-3	Bottle L 100 mL Volume R1 0 mL R2L 0 mL	M 0 mL S 0 mL M 20 mL S								
Rerun (High) 4 uL Dilution <input checked="" type="radio"/> Disable <input type="radio"/> Enable	Rerun (Low) 16 uL	Instrument Factor a 1	Stirring Speed R1 Mid	R2Mid Mid	Abs Range SD								
Save	Cancel	Standard/Ca No.	Expiry										
Calibration (F7)	<input checked="" type="checkbox"/> Check Blank measurement <input type="checkbox"/> Disable reagent blank <input type="radio"/> Enable S1 blank (Factor or linear) <input type="radio"/> Enable reagent blank <input type="radio"/> Enable reagent blank for S1 (linear)	S1: WATER * 0.0 * S2: TruCal * * *	<table border="1"> <thead> <tr> <th>No</th> <th>Value</th> <th>Expiry</th> <th>Abs</th> </tr> </thead> <tbody> <tr> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>	No	Value	Expiry	Abs						
No	Value	Expiry	Abs										
Linearity parameter	□□□□	TruLab N	* * * *	* *	* *								
Factor=	□□□□	TruLab P	* * * *	* *	* *								

NOTE:

These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing the tests.