

## APPLICATION PROCEDURE

REF

Cont.

908510	4 x 20 ml	3 x 20 ml	Reagent 1
	(78 ml)	1 x 18 ml	Reagent 2
908520	4 x 10 ml	3 x 10 ml	Reagent 1
	(39 ml)	1 x 9 ml	Reagent 2

Additionally offered:

909560SV	1 x 1 ml	Homocysteine Calibrator 0 $\mu$ M
909561SV	1 x 1 ml	Homocysteine Cal./Ctrl. approx. 10 $\mu$ M
909562SV	1 x 1 ml	Homocysteine Cal./Ctrl. approx. 33 $\mu$ M
905620	4 x 1 ml	Homocysteine Control Set (4 levels)
908550	5 x 1 ml	Homocysteine Calibrator Set (5 levels)

### BT 2000/3000 Targa

Double reagent application

Test Type	Fixed Time With Start
Serum Starter	(Not Active)
Filters (A/B)	340 / 700
Units	$\mu$ mol/l
Test Method	Multi Point
Test Methodology	UV
Number Of Needle Washes	2 / 2
Delay Time (Sec)	0
Reagent Inc. (Sec)	180 / 150
Reading Time (Sec)	360
Reaction Limit (mABS)	700
Max ABS Delta (mABS)	1500
Reagent Limit (mABS)	1000
Reagents A/B ( $\mu$ l)	240 / 65
Reaction Direction	Decreasing
Reagent Diluent	1:1
Initial ABS (mABS)	2000
Curve Acceptance	100%
Cross Contamination Map	(Not Active)
Automatic Profile	(Not Active)
Rerun Reagent BLK H: M	00:00
Dynamic Blank	(Not Active)
Pathological Repetition	(Not Active)

### SERUM PARAMETERS

Test's Name	Homocysteine
Sample Volume ( $\mu$ l)	13
Dilution Ratio	1:2
Test Limit	60.0
Min. Max. M.	3.0 /27.0
Min. Max. F.	3.0 /27.0
Min. Max. B.	3.0 /27.0

### URINE PARAMETERS

Test's Name	
Urine Volume ( $\mu$ l)	
Test Limit	
Min. Max. M.	---/---
Min. Max. F.	---/---
Min. Max. B.	---/---
Auto Dilution	Not Active
Multi Factor	1
Units Factors	1

**SAMPLE:** Fresh serum, heparin or EDTA plasma (Discard hemolysed, turbid or severely lipemic specimens.

Centrifuge blood samples immediately after collection!

The Homocysteine in the sample is stable at room temperature 4 days, at 0-8°C several weeks, at -20°C > 1year.

### REFERENCE VALUES

**Serum:** Men 6.0 – 15.0  $\mu$ mol/l, Women 3.0 – 12.0  $\mu$ mol/l

Reference values should be used as a guide only. It is recommended that each laboratory should establish its own reference ranges for patient population. The analytical results should be evaluated with other information coming from patient's clinical history.

### CALCULATION:

All values of the determined samples are referred to the generated calibration curve.

**CALIBRATION:** Use kit REF 620 to create the standard calibration curve. The Standard value is lot dependent. The standard value is indicated in the insert of kit REF 620. The calibrators must be placed on the sample plate observing progressive increase of concentrations.

### Standard Calibration Curve:

STD # 1	STD # 2	STD # 3	STD # 4	STD # 5
Level Zero	Level 1	Level 2	Level 3	Level 4

### REAGENT STORAGE AND STABILITY

The data is printed in "Manual Methodology" insert in the kit

### PREPARATION OF REAGENTS

The reagents are liquid and ready to use.

### QUALITY CONTROLS:

Suggested Calibrators and Control Sample:  
Calibrator REF 620., Controls are included in kit

Application data stated in this sheet should be considered as indicative even if coming from routine tests made with Biotechnica's automatic analyzers. They are not validated in

accordance with the directive 98/79-CE. Each laboratory is responsible for confirming the validity of indicated data through a Quality Control procedure.

### BT 1000//BT Plus series

### BT 2500 / BT 3500

Double reagent application

Code for BAR_CODE	Inactive	Inactive
Test Methodology	UV	
Method	Fixed Time	Fixed Time
Kind Of Process	Multi Point	Multi Point
Filters	340 / 700	340 / 700
Reaction Direction	Decreasing	Decreasing
Reagent #1 ( $\mu$ l)	240	200
Reagent #2 ( $\mu$ l)	65	54
Sample Starter	Inactive	Inactive
Delay Time (sec)	0	0
Incubation Time (sec)	180 / 150	180 / 150
Reading Time (sec)	360	360
Unit Serum	$\mu$ mol/l	$\mu$ mol/l
Unit Urine		
Number Of Needle Washes	2 / 2	2 / 2
Number Of Cuvette Washes	2	2
Dynamic Blank	Inactive	Inactive
Reagent Blank	Every Day	Every Day
Reagent Limit (mABS)	1000	1000
Curves Acceptance (%)	100	100

### SERUM PARAMETERS

Name	Homocysteine	
Sample ( $\mu$ l)	13	10.8
Pre-Dilution	1	1
Dilution		
Factor	2	2
Limit Test (Conc.)	60.0	60.0
Initial ABS (mABS)	2000	2000
Final ABS (mABS)	700	700
Max ABS Delta (mABS)	1500	1500
Check Prozona (mABS)	Inactive	Inactive
Instrumental Factor	1	1
Shift	0.000	0.000
Re-run Hyperactive	Active	Active
Re-run Pathological	Inactive	Inactive
RE-run out of curve "Above"	Active	Active
Re-run out of curve "Below"	Inactive	Inactive
Normal Range		
Men	3.0 /27.0	
Women	3.0 /27.0	
Children	3.0 /27.0	

### URINE PARAMETERS

Name		
Sample ( $\mu$ l)	3	3
Pre-Dilution	1	1
Diluzione		
Factor	1	1
Limit Test (Conc.)		
Initial ABS (mABS)		
Final ABS (mABS)		
Max ABS Delta (mABS)		
Check Prozona (mABS)	Inactive	Inactive
Instrumental Factor	1	1
Shift	0.000	0.000
Re-run Hyperactive	Inactive	Inactive
Re-run Pathological	Inactive	Inactive
RE-run out of curve "Above"	Inactive	Inactive
Re-run out of curve "Below"	Inactive	Inactive
Normal Range		
Men		
Women		
Children		

### INTERFERENCES

The data is printed in "Manual Methodology" insert in the kit.

### LINEARITY

The data is printed in "Manual Methodology" insert in the kit. The declared data can be obtained with the instrument under optimal conditions of use.

**ANALYTICAL PERFORMANCE:** The data is printed in "Manual Methodology" insert in the kit. The declared data can be obtained with the instrument under optimal conditions of use.

### NOTE

- Avoid reagent exposure to heating, direct sunlight, contamination and evaporation. In case of complaint or quality control request, indicate the lot number of the package or the lot number of the individual components.

**WASTE MANAGEMENT:** Dispose of in accordance with local official regulations.

	Consult accompanying documents		Consult instructions for use		Manufactured by DIALAB GmbH
	Biological risks		Do not dispose of in environment		Rev01. 09.04.2010