

HITACHI 917

CHOLESTEROL LDL, DIRECT ENZYMATIC SELECTIVE PROTECTION

Diagnostic reagent for quantitative in vitro determination of low density lipoprotein cholesterol (LDL-C) in human serum or plasma on photometric systems

REF	Cont.				
F05365	5 x 50 ml	4 x 50 ml	Reagent 1	1250 tests/kit	
		1 x 50 ml	Reagent 2		
F05367	5 x 10 ml	4 x 10 ml	Reagent 1	250 tests/kit	
		1 x 10 ml	Reagent 2		
F0443917*	5 x 62,5 ml	4 x 62,5 ml	Reagent 1	1562 tests/kit	
		1 x 62,5 ml	Reagent 2		
F0417917*	5 x 20 ml	4 x 20 ml	Reagent 1	500 tests/kit	
		1 x 20 ml	Reagent 2		

* Hitachi System Pack

Additionally offered:

F03711SV	1 x 1 ml	LDL Cholesterol Calibrator
D99486	3 x 3 ml	Lipid Control normal Diacon Lipids

1. Reagent preparation:

The reagents are ready to use.

2. Instrument settings:

Analysis		LDL				Ser/PI			
Test / Type		2	A	10	A	15	33	0	0
Assay / Time / Point		700	A	600	A				
Wave (2 nd /Primary)									
S.Vol (Normal)		2	0	0					
S.Vol (Decrease)		2	0	0					
S.Vol (Increase)		5	0	0					
Diluent		00951	99						
Reagent (R1) T1		160	0	059	0				R1
Reagent (R2) T2		0	0	059	0				R2
Reagent (R3) T3		40	0	059	0				R3
Reagent (R4) T4		0	0	059	0				
Abs. Limit		12000	Increase		A				
Prozone Limit		32000	0	Lower	A				
Cell Detergent		Detergent 1		A					

Calibration		Linear		A		A	
Calibration type		2	Span Point	2			
Point							
Weight		0					
Autocalibration							
Time Out				Change Over			
Blank				Blank		A	
Span				Blank		A	
2Point							
Full							
SD Limit		0.1					
Duplicate limit		10	%	200		Abs	
Sensitivity limit		-99999		999999			
S1 Abs limit		-32000		32000			

APPLICATION PROCEDURE



Range			
Application Code	059		Unit mg/dl A
Report Name	LDL C		
Data Mode	On Board	A	
Control Interval	1000		
Instrument Factor (Y=aX+b)	a= 1.0	b=	0.0
Technical Limit	0	400	
Repeat Limit	0	400	
Expected Value			
Qualitative			
(Male)	Y	A	
	Y	A	
			0 130
(Female)	Y	A	
	Y	A	
(Default)	Male	A	Range3 A
			(1) 0
			(2) 0
			(3) 0
			(4) 0
			(5) 0
			(6) 0

Others						
<Standard>	(1)	(2)	(3)	(4)	(5)	(6)
Calib. Code	501	#				
Concentration	0.0	*				
Position						
Sample Volume	2	2				
Diluent S.Vol	0	0				
Diluent Volume	0	0				

- #) Data entry by the user
- *) Enter calibration or standard value

NOTE: These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing the tests.