

## ALCYON-FALCOR

### $\alpha$ -AMYLASE

mod. IFCC

Liquid

2 Reagents

REF

Cont.

<b>D94570</b>	<b>5 x 100 ml</b>	4 x 100 ml	Reagent 1	<b>2500 tests</b>
		1 x 100 ml	Reagent 2	
<b>D94571</b>	<b>5 x 50 ml</b>	4 x 50 ml	Reagent 1	<b>1250 tests</b>
		1 x 50 ml	Reagent 2	
<b>D96569</b>	<b>5 x 10 ml</b>	4 x 10 ml	Reagent 1	<b>250 tests</b>
		1 x 10 ml	Reagent 2	

D98485	5 x 3 ml	Calibrator	Diacal Auto
D98481	12 x 5 ml	Control normal	Diacon N
D98482	12 x 5 ml	Control abnormal	Diacon P

#### 1. Reagent preparation:

The reagents are ready to use. For Sample Start mix 4 parts of Reagent 1 with 1 part of Reagent 2.

#### 2. Instrument settings:

Type	: Kin 0	Primary Wavelength	: 405
Unit	: U/l	High Wavelength	: -
Number of Decimals	: 1	Low Wavelength	: -
Diluent	: 99 Diluan		:
Pre-washing	: No	High Limit	: 0.07
First reagent Vol. ( $\mu$ l)	: 200	Depletion limit	: 0.2400
Second reagent Vol. ( $\mu$ l)	: 0	C Compensation factor	: 0.0650
Sample Volume ( $\mu$ l)	: 4.0		
Lower normal value	: 0	Slope	: 1.00
Higher normal value	: 100	Intercept	: 0.0
Dilution factor	: 4		
Systematic dilution	: 1		
Dilution for Urine	: 1		
Low reagent blank	: 0.0000	Reaction delay	: 02:00
High reagent blank	: 0.3000		
Maximum Deviation	: 0.0300		

#) Data entry by the user

\*) Calculated by the analyzer

#### NOTE:

These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing the tests.