



FURUNO 180 / DAYTONA

GLUCOSE

GOD - PAP

Liquid, Single Reagent

REF

Cont.

D95218B	1 x 1000 ml	Single Reagent	5000 tests/kit
D08220	4 x 250 ml	Single Reagent	5000 tests/kit
D00221	5 x 100 ml	Single Reagent	2500 tests/kit
D98485	5 x 3 ml	Calibrator	Diacal Auto
D98481	12 x 5 ml	Control normal	Diacon N
D98482	12 x 5 ml	Control abnormal	Diacon P

1. Reagent preparation:

The reagent is ready to use.

2. Instrument settings:

RunMnr(F5)	ChmPrm(F6)	Calibration (F7)	QC (F8)	System Prms	Maintenance
Chemistry Parameters Method: 14 Name: GLUC		R1 Reagent Name: GLUC R2: <input type="radio"/> Disable <input checked="" type="radio"/> Enable		Volume: 200 uL	
Unit: mmol/L Assay Type: End		Reagent Name: _____ Wash: <input type="radio"/> Disable <input checked="" type="radio"/> Enable		Volume: _____ uL Reagent Name: _____ Reagent Type: Wash	
Measuring Points: 1 <input checked="" type="radio"/> Disable	Start: - End: -	Dilution: <input type="radio"/> Disable <input checked="" type="radio"/> Enable	Start: 25 End: 26	Decimal Points: 1 Normal Range: 3.9 - 6.4	Reagent Name: _____
Wave Length: Prim. 510 Sec. <input type="checkbox"/>				Extension: _____	
Sampling Volume: 2 uL Dilution: <input checked="" type="radio"/> Disable <input type="radio"/> Enable Diluent: 0 uL		Technical Range (Conc.): 0 - 22.2 (mAbs/10)-50000 - 50000 RPT Wash (R1): <input checked="" type="radio"/> SysWater <input type="radio"/> Wash Sol-3 (R2): <input checked="" type="radio"/> SysWater <input type="radio"/> Wash Sol-3			
Rerun (High): 2 uL Dilution: <input checked="" type="radio"/> Disable <input type="radio"/> Enable Diluent: 0 uL		Bottle Volume L: 100 mL R1: 0 mL R2L: 0 mL		Instrument Factor a: 1 b: 0	
Rerun (Low): 4 uL				Stirring Speed R1: Mid Standard/Cal No. value. Expiry	R2 Mid: _____ Abs Range: SD
Calibration (F7) Checks (F10)	Blank measurement <input type="checkbox"/> Disable reagent blank and S1 blank <input checked="" type="radio"/> Enable S1 blank (Factor or Linear) <input type="radio"/> Enable reagent blank <input type="radio"/> Enable reagent blank for S1 (Linear)	S1: WATER * 0.0 * S2: TruCal U * * *			
Linearity parameter=	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	TruLab N * * *			
Factor=	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	TruLab P * * *			

NOTE:

These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing the tests.