

TARGA BT 3000

GLUCOSE GOD-PAP

Liquid
Single Reagent

REF

Cont.

D95218B	1 x 1000 ml	Single Reagent	3333 tests 300 µl/test
D00221	5 x 100 ml	Single Reagent	1666 tests 300 µl/test

D95223	1 x 3 ml	Glucose Standard	
D98485	5 x 3 ml	Calibrator	Diacal Auto
D98481	12 x 5 ml	Control normal	Diacon N
D98482	12 x 5 ml	Control abnormal	Diacon P

1. Reagent preparation:

The reagent is ready for use.

2. Instrument settings:

Test Typ	Endpoint Single	Serum Parameters	
Serum Starter	(Inactive)	Tests Name	GLUCOSE
Filters (A/B)	510/700	Sample Volume	(µl) 3
Units	mg/dl	Dilution Ratio	1:31
Test Method	With Factor	Min. Max. M.	75/115
Test Methodology	TRINDER	Min. Max. F.	0.000/0.000
Number of washes	1	Min. Max. B.	0.000/0.000
Delay Time	(Sec) 0	Urine Parameters	
Inc. Time	(Sec) 600	Tests Name	
Reading Time	(Sec) 10	Urine Volume	(µl) 3
Test Limit	(Conc) 400	Min. Max. M.	0.000 / 0.000
Reactio Limit	(mABS) 0.000	Min. Max. F.	0.000 / 0.000
Max ABS Delta	(mABS) 0.000	Min. Max. B.	0.000 / 0.000
Reagent mAbs Limit	(mABS) 200	Auto Dilution	1:3
Reagents A/B	(µl) 300	Multi Factor	1.00
Reaction Direction	Increasing		
Reagent Dilution	1:1		
Initial ABS	(mABS) 0.000		
Curve Acceptance	100%		
Automatic profile	(Inactive)		
Rerun Test Rgt Blk	H:M 00:00		
Pathological Repetition	(Inactive)		

NOTE:

These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing the tests.