

APPLICATION PROCEDURE

MaxMat PL

GPT (ALT, ALAT)

MODIFIED IFCC

Liquid
2 Reagents

REF

Cont.

D94620	5 x 100 ml	4 x 100 ml 1 x 100 ml	Reagent 1 Reagent 2	3125 tests 160 µl/test
D98624	5 x 50 ml	4 x 50 ml 1 x 50 ml	Reagent 1 Reagent 2	1562 tests 160 µl/test

D98485	5 x 3 ml	Calibrator	Diacal Auto
D98481	12 x 5 ml	Control normal	Diacon N
D98482	12 x 5 ml	Control abnormal	Diacon P

1. Reagent preparation:

The Reagents are ready to use. For Sample Start mix thoroughly 4 parts of Reagent 1 with 1 part of Reagent

2. Instrument settings:

Sample Start:

Temperature: 37 °C

METHOD PARAMETERS					EXPECTED VALUES AND LIMITS	
DETERMINATION TYPE:					EXPECTED VALUES	
START POINT IS THE MEAN OF					MIN. VALUE	0
END POINT IS THE MEAN OF					MAX VALUE	41
REACTION WAY:					METHODIC DYNAMIC RANGE	
CHROMATISM:					LOW VALUE	4
MAIN:					HIGH VALUE	349
ASSOCIATED 1:					REAGENT ABSORBANCE LIMIT	NO
ASSOCIATED 2:					SUBSTRATE DEPLETION LIMIT	YES
SHOW WAVELENGTHS						CYCLE 4
POLYCHROMATIC CORRECTION:					BICHROMATIC ALARM	NO
CORRECTION FACTOR :					CORRELLATION FACTOR	
MAIN UNIT:					A	0
SECOND UNIT					B	1
RESULT FORMAT:					AUTOMATIC RERUN	NO (OR #)
PRE-DILUTION:					POST CONCENTRATION	NO (OR #)
RATIO:					POST DILUTION	YES
DILUENT:					RATIO	1/5
RELATIVE CORRECTION FACTOR:					DILUENT	Chlorure de sodium 9 g/l
CONNECTION DILUENT						
METHOD STEPS					CALIBRATION	
CYCLE	ACTION	VOLUME (µL)	REAGENT	FLUSH VOL	CALCULATION CURVE FITTING	LINEAR
1	REAGENT	160	ALT WORKING R.	15	CALIBRATION FACTOR	-
2	SAMPLE	20		15	STANDARD	DIACAL A (OR#)
MEASURE: NORMAL; HIGH SPEED MIXING					STANDARD DUPLICATE	YES
3-16	NOTHING				MANUAL CALIBRATION	NO
6	START MEASURE CYCLE				AUTOMATIC VALIDATION	NO
15	END OF MEASUREMENT				NO DILUTION	STANDARD 1/1
16	END OF METHOD				RELATIVE CORRECTION FACT.	1

#) Data entry by the user

NOTE:

These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing the tests.