

ABX PENTRA 400

GPT (ALT, ALAT) MODIFIED IFCC

Liquid
2 Reagents

REF Cont.

D94620	5 x 100 ml	4 x 100 ml	Reagent 1	2500 tests
		1 x 100 ml	Reagent 2	
D98624	5 x 50 ml	4 x 50 ml	Reagent 1	1250 tests
		1 x 50 ml	Reagent 2	

D98485	5 x 3 ml	Calibrator	Diacal Auto
D98481	12 x 5 ml	Control normal	Diacon N
D98482	12 x 5 ml	Control abnormal	Diacon P

1. Reagent preparation:
The Reagents are ready to use.

2. Instrument settings:

GENERAL PARAMETERS			
Test Name	Code	Channel	Release
<input type="text" value="ALT"/>	<input type="text" value="ALT"/>	<input type="text" value="#"/>	<input type="text" value="App. Version"/>
CHARACTERISTICS		<input type="checkbox"/> Pre-Dilution <input checked="" type="checkbox"/> Enable	<input type="text" value="#"/> Modified on
Sample Type	<input type="text" value="Serum/Plasma"/>	Diluent Name	<input type="text" value="U/l"/> Unit
Number of Reagents	<input type="text" value="Reagent 2"/>	Factor	<input type="text" value="0"/> Decimal position
REAGENT		Incubation time (cycles)	<input type="checkbox"/> Manual Patient Validation
Reagent short name	Reagent number		
<input type="text" value="ALTR"/>	<input type="text" value="#"/>		
On Board Stability	<input type="text" value="28"/>	Linearity	Correlation
#	Cassette	Low limit	Slope
		<input type="text" value="4.00"/>	<input type="text" value="1.0000"/>
		High limit	Intercept
		<input type="text" value="600.00"/>	<input type="text" value="0.0000"/>
Delta Check			
Automatic Rerun		Delta Check Validity	Absolute Variation
<input checked="" type="checkbox"/> Post Dilution	Dilution Factor <input type="text" value="3.0"/>	<input type="text"/>	Relative Variation
<input type="checkbox"/> Post Concentration	Concentration Factor <input type="text"/>	Reference Range	
		<input checked="" type="checkbox"/> Low Check	Man/Default
		<input type="text" value="0"/>	Woman
		<input checked="" type="checkbox"/> High Check	Child
		<input type="text" value="41.00"/>	<input type="text"/>
		<input type="text"/>	<input type="text"/>

CALIBRATION PARAMETERS			
Test Name	Code	Channel	Checks
<input type="text" value="ALT"/>	<input type="text" value="ALT"/>	<input type="text" value="#"/>	<input checked="" type="checkbox"/> Reagent Limit Absorbance Check
Pre-Dilution			Reagent Range Low <input type="text" value="-0.2000"/>
Type	Calibrator	Diluent	Reagent Range High <input type="text" value="1.6000"/>
Factor 1	Factor 2	Factor 3	<input checked="" type="checkbox"/> Reagent Blank Limit Absorbance Check
<input type="text"/>	<input type="text"/>	<input type="text"/>	Blank Range – Low Limit <input type="text" value="-0.0400"/>
Factor 5	Factor 6	Factor 7	Blank Range – High Limit <input type="text" value="0.0100"/>
<input type="text"/>	<input type="text"/>	<input type="text"/>	
Factor 8	<input type="text"/>	<input type="text"/>	
<input type="text"/>			



APPLICATION PROCEDURE

Calibration Calibration Mode	Slope Avg	Validity <input checked="" type="radio"/> On Request <input type="radio"/> Time Validity Interval Time Unit # Days	Control Required Control Used _____
Level	1		
Calibration Factor	*		
Runs	2		
<input checked="" type="checkbox"/> Deviation %	5.0		
Calibrator Used			

ANALYSIS PARAMETERS						
Test Name	Code	Channel				
ALT	ALT	#				
<input type="checkbox"/> Cleaner	Wavelength (nm)		Blank			
Cleaner Solution	H2O	Primary Wavelength	340	<input checked="" type="checkbox"/> Reagent Blank	H2O	Mixing Speed 80
<input type="radio"/> Before		Secondary Wavelength	420			
<input type="radio"/> After						
Analysis Sequence						
Cycle	Reagent Needle	Volume (µl)	Sample Needle	Volume (µl)	H2O Vol (µl)	
1	R1	160.0	SAMPLE	20.0	10.0	
15			R2	40.0	10.0	

CALCULATION PARAMETERS						
Test Name	Code	Channel				
ALT	ALT	#				
Correlation Factor	Reaction Direction		<input checked="" type="checkbox"/> Sample Limit Check			
Slope	1.0000	<input checked="" type="checkbox"/> Reaction Direction Check	Decrease	Sample limit (Δ O.D.)	0.1200	
Intercept	0.0000	Reaction Direction		Sample limit cycle	2	
Definition						
Calculation Type	Kinetic		First Reading		Last Reading	
<input checked="" type="checkbox"/> Reaction limit Check	Reaction limit Absorbance		Cycle	18	Cycle	30
	0.0800					
Cycle	16					

UNIT PARAMETERS	
Test Name	ALT
Code	ALT
Channel	#
Unit	Conversion Factor
U/I	1.0000
Modified On	

NOTE:
These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing the tests.